

Datasheet

FUBP3 mouse monoclonal antibody (hybridoma)

Catalog Number: H00008939-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant FUBP3.

Immunogen: FUBP3 (AAH07874.1, 1 a.a. ~ 261 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MPPIGNQLGALVHQRTVITEEFKVPDKMVGFIIGRGGE
QISRIQAESGCKIQIASSESGIPERPCVLTGTPESIEQAK
RLLGQIVDRCRNGPGFHNDIDSNSTIQEILIPASKVGLVI
GRGGETIKQLQERTGVKMVMIQDGPLTGADKPLRIT
GDAFKVQQAREMVLEIIREKDQADFRGVRGDFNSRM
GGGSIEVSVPRFAVGIVIGRNGEMIKKIQNDAGVRIQFK
PDDGISPEYYRQQVAFYGGTLGQAQAHSEQ

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Entrez GeneID: 8939

Gene Symbol: FUBP3

Gene Alias: FBP3