

Datasheet

SMC2L1 polyclonal antibody (A01)

Catalog Number: H00010592-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant SMC2L1.

Immunogen: SMC2L1 (NP_006435, 521 a.a. ~ 630 a.a) partial recombinant protein with GST tag.

Sequence:

CVKGLVASLISVKDTSATTALELVAGERLYNVVVDTEV
TGKKLLERGELKRRYTIPLNKISARCIAPETLRVAQNLV
GPDNVHVALSLVEYKPELQKAMEFVFGTTFVC

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10592

Gene Symbol: SMC2

Gene Alias: CAP-E, CAPE, FLJ10093, SMC2L1, hCAP-E